ADJUSTABLE SLEEP & CARE SOLUTIONS

Supporting mobility, positioning, comfort and independence—without compromising the feeling of home.



DESIGNER DECOR

ENHANCES COMFORT

INDEPENDENT LIVING

ASSISTS CAREGIVERS

Designed for home care, rehabilitation, aged care and everyday living.

THERAPEDIC® WORLDWIDE



DESIGNED AROUND THE PERSON

Care-focused functionality. Residential comfort.

The right sleep environment can play an important role in supporting comfort, mobility and independence at home. Therapedic adjustable beds and specialised sleep systems are designed to assist with everyday activities including bed entry and exit, repositioning, supported sitting, resting and sleeping.
With a choice of adjustable functions, bed configurations, mattress feels and support options, Therapedic makes it easier to select a solution around the individual—rather than expecting the individual to adapt to the bed.
Our approach combines practical care functionality with contemporary furniture styling, helping create a supportive sleep environment that still feels like home.

FOR HEALTH PROFESSIONALS

Therapedic provides detailed product specifications, dimensions, mattress options and configuration information to assist Occupational Therapists and other healthcare professionals when considering an adjustable sleep solution for an individual client.

OUR COMMITMENT

To provide dependable sleep and mobility solutions that support independence, everyday function, comfort and quality of life.



CONTENTS

PARTNER BED OPTIONS P.1
BED SIZE & DIMENSIONS P.2
390 ADJUSTABLE BED & POSITIONS P.3-P.4
460 BEDS & FUNCTIONS P.5-P.6
THERAPEDIC BED & POSITIONS P.7-P.10
THERAFLEX T-SERIES MATTRESS P.11-P.12
T15 MATTRESS SPECIFICATIONS P.13-P.14
T20 MATTRESS SPECIFICATIONS P.15-P.16
T25 MATTRESS SPECIFICATIONS P.17-P.18
T30 MATTRESS SPECIFICATIONS P.19-P.20
PROFILE SPRING P.21-P.22
PROFILE THERAPY COMFORT MATTRESS . P.23-P.24
PROFILE THERAPY LUXURY MATTRESS P.25-P.26

PARTNER BED OPTIONS

For couples, a Partner Bed configuration can allow an individual to access adjustable functionality while maintaining the appearance and familiarity of a shared bed. This can be particularly valuable when the support needs of one person change while their partner's do not.

CARE THAT STILL FEELS LIKE HOME

A change in mobility or support needs doesn't necessarily mean giving up the comfort, appearance or familiarity of a conventional bedroom.
Therapedic adjustable beds combine powered positioning and practical care functionality with contemporary residential design, allowing the bed to integrate naturally into the home.



CHANGING MOBILITY NEEDS

Bed configurations that can be selected around functional requirements.



COUPLES

Options designed to preserve a shared bedroom environment.




THE HOME

Healthcare functionality without an institutional appearance.

BED SIZE & DIMENSIONS

	Split Queen	Long Single	King Single	Double	Queen	Partner Base (all sizes)
	77cmx203cm	92cmx203cm	107cmx203cm	138cmx188cm	153cmx203cm	
Theracloud		✔	✔	✔	✔	✔
Theralift Elite	✔	✔	✔	✔	✔	✔
Theramotion Advance		✔	✔	✔	✔	✔
SBAD0300460	✔	✔	✔	✔		
SBAD0300390		97cmx208cm Inc.surround	113cmx208cm Inc.surround	143cmx193cm Inc.surround		
		✔	✔	✔		




Clinical consideration: Bed size should be considered alongside the individual's mobility, transfer technique, positioning requirements, available bedroom space and any equipment required around the bed.




390 HI-LO ADJUSTABLE BED & POSITIONS

The Therapedic 390 combines powered head and leg positioning with Hi-Lo height adjustment, providing greater flexibility for sleeping, resting, transfers and day-to-day care.




UPPER BODY ELEVATION

Anti-Snore preset – raises the head and upper body to provide an elevated resting or sleeping position.




SUPPORTED REST

Zero Gravity preset – coordinates upper and lower body elevation to create a relaxed, supported resting position and redistribute load.




SUPPORTED SITTING

Raises the head and upper back for activities such as reading, watching television, communicating or relaxing in bed.




FLAT POSITION

Returns the mattress platform to a conventional flat sleep position.



LOWER-LIMB ELEVATION

Raises the lower section of the bed independently for personalised positioning and comfort.



HI-LO HEIGHT ADJUSTMENT

Raises or lowers the complete sleep surface to suit transfers, bed access and care activities.

OT / Functional consideration

Consider where variable bed height may assist with the client's transfer strategy or caregiver access.

IMPORTANT

Guard rails and bedside supports should only be selected following consideration of the user's mobility, transfer method, mattress configuration and individual risk factors.

Weight Capacity: 220KG

Bed Height:

(Flat not inc. mattress)

Max Height : 62cm
Minimum Height: 33cm

2 Years Warranty Remote | 6 Years Warranty Parts | 12 Years Warranty Frame

BEDS & FUNCTIONS

Smart Features for Better Living

Engineered for effortless daily living, adjustable beds combine intuitive controls, smooth Hi-Lo height adjustment, and thoughtful safety features like guide rails and mattress retainers for unmatched convenience, care, and comfort.



Wired Remote

Intuitive tactile buttons designed for effortless, fail-safe operation—keeping core adjustments within arm's reach without battery concerns.



Lockable Wheels

High-load, smooth-rolling castors allow seamless bed repositioning for cleaning or care, tucking neatly beneath the surround with rock-solid locking.



Easy Clean Metal Grid

Powder-coated steel grid structure maximizes airflow beneath the mattress while enabling quick, thorough sanitisation for daily care.



Guard Rail

Engineered for safety and transfer assistance, blending clean architectural lines with dependable support to prevent accidental falls.



Hi-Lo Adjustable Base

Smooth, powered vertical travel optimizes ergonomics for caregivers and aligns perfectly with individual transfer heights.



Removable Side Fabric

High-grade upholstery attaches cleanly to the surround, allowing quick removal for laundering without compromising the bed's tailored aesthetic.




End Bar - Stop Mattress Moving

Solid stainless-steel retention bracket holds the mattress precisely in position throughout articulate head, leg, and height movements.




460 HI-LO ADJUSTABLE BED & POSITIONS

The Therapedic 460 combines powered positioning and Hi-Lo height adjustment with a fully upholstered surround designed to sit comfortably within a contemporary bedroom.




UPPER BODY ELEVATION

Anti-Snore preset – raises the head and upper body to provide an elevated resting or sleeping position.




SUPPORTED REST

Zero Gravity preset – coordinates upper and lower body elevation to create a relaxed, supported resting position and redistribute load.




SUPPORTED SITTING

Raises the head and upper back for activities such as reading, watching television, communicating or relaxing in bed.




FLAT POSITION

Returns the mattress platform to a conventional flat sleep position.



LOWER-LIMB ELEVATION

Raises the lower section of the bed independently for personalised positioning and comfort.



HI-LO HEIGHT ADJUSTMENT

Raises or lowers the complete sleep surface to suit transfers, bed access and care activities.

OT / Functional consideration

Consider where Hi-Lo functionality is required but maintaining a residential bedroom appearance is also important.

Weight Capacity: 220KG

Bed Height:

(Flat not inc. mattress)

Max Height : 67cm
Minimum Height: 38cm

BEDS & FUNCTIONS

Smart Features for Better Living

Designed to blend luxury and function—featuring intuitive wireless positioning, relaxing built-in massage, non-slip fabric, and elegant foam cushioning for a superior sleeping and relaxation experience.




Wireless Remote

Ergonomically contoured with responsive tactile feedback, featuring one-touch lifestyle presets for seamless, cable-free comfort control.



Lockable Wheels

High-load, smooth-rolling castors allow seamless bed repositioning for cleaning or care, tucking neatly beneath the surround with rock-solid locking.



Non-Slip Fabric

High-friction technical textile across the articulation platform secures the mattress seamlessly, preventing shifting even at maximum elevation angles.




Fully Upholstered Foam Surround

Premium tailored fabric over dense, protective foam cushioning softens impact, transforming clinical functionality into contemporary bedroom elegance.



Hi-Lo Adjustable Base

Smooth, powered vertical travel optimizes ergonomics for caregivers and aligns perfectly with individual transfer heights.



Dual-Zone / Built-In Massage

Integrated whisper-quiet massage motors deliver soothing vibrational therapy across targeted zones to ease muscle tension and promote restorative circulation.



End Bar - Stop Mattress Moving

Solid stainless-steel retention bracket holds the mattress precisely in position throughout articulate head, leg, and height movements.




THERAPEDIC ADJUSTABLE RANGE
THERACLOUD, THERALIFT ELITE & THERAMOTION ADVANCE

Therapedic adjustable bases bring powered positioning into a contemporary residential bed.




UPPER BODY ELEVATION

Anti-Snore preset – raises the head and upper body to provide an elevated resting or sleeping position.




SUPPORTED REST

Zero Gravity preset – coordinates upper and lower body elevation to create a relaxed, supported resting position and redistribute load.




SUPPORTED SITTING

Raises the head and upper back for activities such as reading, watching television, communicating or relaxing in bed.




HEAD TILT POSITION (THERAMOTION ADVANCE ONLY)

Adjusts the head section forward and upward independently from the backrest. This targeted movement fills the gap beneath the neck and provides optimal cervical support, allowing you to comfortably read, watch TV.



LOWER-LIMB ELEVATION

Raises the lower section of the bed independently for personalised positioning and comfort.



FLAT POSITION

Returns the mattress platform to a conventional flat sleep position.

OT / Functional consideration

Consider whether independent head and leg positioning may improve the user's ability to self-position for rest and activities in bed.

Weight Capacity: 220KG

Bed Height:

(Flat not inc. mattress)

Max Height : 40cm
Minimum Height: 12cm

BEDS & FUNCTIONS

Smart Features for Better Living

The innovative GEO Curve design features two closely positioned articulation points, creating a gentler curve compared to traditional adjustable beds with a single hinge. This allows the mattress to flex more naturally, reducing material stress while delivering improved body conformity, pressure relief, and overall comfort. With modern smart tech—including targeted Head Tilt, relaxing massage, ambient LED lighting, and USB charging—for personalized wellness at the touch of a button.



Adjustable Bed Legs

Stackable and detachable, this heavy-duty steel leg lets you easily customize your bed height while ensuring solid stability for adjustable bases.



Wireless Remote

Ergonomically contoured with responsive tactile feedback, featuring one-touch lifestyle presets for seamless, cable-free comfort control.



LED Underbed Light

Soft under-bed ambient illumination casts a warm, glare-free halo across the floor, ensuring safe nocturnal navigation without disturbing a resting partner.



USB Charging Ports

Discretely positioned high-speed charging ports keep smartphones, tablets, and bedside essentials powered directly from the comfort of your bed.



Head Tilt (Theramotion Advance Only)

Independent adjustment of the upper head section provides additional support for your head and neck, helping you find a more comfortable position for reading, watching TV, or relaxing.




Massage (Theralift Elite & Theramotion Advance Only)

Whisper-quiet therapeutic vibration motors deliver pulsating wave massage across multiple intensity levels, alleviating muscle fatigue and encouraging restorative sleep

GEO CURVE (AUTO-REGRESSION)

(CONTOURED COMFORT, THROUGH EVERY MOVEMENT)


A unique feature across the range, this subtle contouring supports your spine's natural curve when elevating. Unlike flat or rigid frames, Geo-Curve helps maintain neutral spinal alignment, easing back tension and improving comfort, especially when paired with zero-gravity or head-raise positions.



CONVENTIONAL ARTICULATION

Single-point articulation

A larger gap can form between the mattress and base, creating a more angular profile as the bed moves.



GEO CURVE ARTICULATION

A gentler, closely supported curve

Multiple, closely positioned articulation points allow the mattress to flex naturally, maintaining closer contact with the base and supporting a smoother profile.

- **German-Designed Okin Motors -** Ultra-quiet and precision-built for long-lasting use.
- **Semi-Fowler's Positioning -** Elevate the upper body 30–45° to relieve snoring, reflux, breathing issues, and aid circulation, a hospitalrecommended.
- **Custom Elevation -** Independent head and foot adjustment to reduce swelling, relieve joint pain, and enhance blood flow.
- **Engineered for Everyday Performance**
- **320kg Weight Capacity -** Commercial-grade strength designed to support more bodies, more safely.
- **12-Year Frame Warranty -** Durability you can count on, backed by extended manufacturer support.

Feature	Theracloud	Theralift Elite	Theramotion Advance
OKIN MOTORS	2	2	3
Head & foot adjustment	✔	✔	✔
Geo Curve	✔	✔	✔
Head Tilt	✘	✘	✔
Underbed LED	✔	✔	✔
Massage	✘	✔	✔
Dual USB charging ports	✔	✔	✔
Bluetooth app	✔	✔	✔
Alarm Function	✔	✔	✔
Commercial Grade Fabric	✔	✔	✔
Foldable Package	✔	✔	✔

COMMERCIAL
GRADE FABRIC

320 KG
WEIGHT CAPACITY

FOLDABLE DESIGN
Folds for easier delivery and setup

2 Years Warranty Remote

6 Years Warranty Parts

12 Years Warranty Frame



THERAFLEX T-SERIES SUPPORT MATTRESSES


The Mattress Matters

An adjustable bed performs best when paired with a mattress designed to articulate with it.

The Theraflex T-Series—T15, T20, T25 and T30—has been engineered to flex and conform with Therapedic adjustable bases while maintaining structural support.

Across the range, layered foam constructions provide different comfort profiles, allowing the mattress to be selected around the user's comfort preference, positioning requirements and support needs.

THERAFLEX T-SERIES MATTRESS







TheraFlex

OT consideration: Encourage adults who have been assessed as being at risk of developing pressure ulcers to change their position frequently and at least every 6 hours. Encourage adults who have been assessed as being at high risk of developing a pressure ulcer to change their position frequently and at least every 4 hours.

T15

Available in all sizes


Theraflex T15

Feel Rating: Firm

User Weight: Max 320kg


Product Highlight:

- Designed for use on adjustable beds
- Flip free design
- Anti slip base layer
- Compatible with waterproof covers



PRESSURE MAP

Powered & Tested By:
XSENSOR® Technology Corporation
Innovators in Pressure Imaging



TheraFlex T Series Mattress

Standard Mattress

Pressure Gauge

Extreme High Medium Low Minimum

The way a support surface distributes body weight is an important consideration when selecting a mattress for users with individual positioning or pressure-management requirements.

Theraflex mattress uses layered premium foams designed to distribute load across the mattress surface while remaining flexible enough to articulate with the adjustable base.

Pressure mapping is a comparative assessment tool. Individual outcomes vary according to body shape, weight, position, mobility and other factors. Mattress selection should form part of an individual assessment where pressure-injury risk is present.

1. Edge Support (65mm Width)

Extra firm edge integrated throughout the mattress to provide complete structural stability and prevent roll-out.

2. Contact Comfort Layer

Soy Visco memory foam, the most natural foam in the market,one of the main ingredients is soy oil from sustainable plantations. It features incredible elasticity and breathability.

3. Pressure Reduction Layer

High density premium foam provides a firm comfort level.

4. Premium Support Layer

Engineered with high-density foam, the bottom support layer offers superior structural integrity and weight distribution. It is designed to withstand frequent patient use while minimizing mattress sagging.

5. Airflow Fabric

Promotes breathability keeping the mattress cooler.

THERAPEDIC® WORLDWIDE

SOY VISCO®

PREMIUM FOAM

EDGE SUPPORT

AUSTRALIAN MADE

10 years Guarantee

TheraFlex

OT consideration: Encourage adults who have been assessed as being at risk of developing pressure ulcers to change their position frequently and at least every 6 hours. Encourage adults who have been assessed as being at high risk of developing a pressure ulcer to change their position frequently and at least every 4 hours.

T20

Available in all sizes


Theraflex T20

Feel Rating:
Available in Firm, Medium & Plush

User Weight:
Max 320kg

Product Highlight:

- Designed for use on adjustable beds
- Flip free design
- Anti slip base layer
- Compatible with waterproof covers



1. Edge Support (65mm Width)

Extra firm edge integrated throughout the mattress to provide complete structural stability and prevent roll-out.

2. Contact Comfort Layer

Soy Visco memory foam, the most natural foam in the market,one of the main ingredients is soy oil from sustainable plantations. It features incredible elasticity and breathability.

3. Premium Support Layer


Engineered with high-density foam, the bottom support layer offers superior structural integrity and weight distribution. It is designed to withstand frequent patient use while minimizing mattress sagging.

4. Airflow Fabric

Promotes breathability keeping the mattress cooler.

PRESSURE MAP

Powered & Tested By:
XSENSOR® Technology Corporation
Innovators in Pressure Imaging



TheraFlex T Series Mattress

Standard Mattress

Pressure Gauge

Extreme | High | Medium | Low | Minimum

The way a support surface distributes body weight is an important consideration when selecting a mattress for users with individual positioning or pressure-management requirements.

Theraflex mattress uses layered premium foams designed to distribute load across the mattress surface while remaining flexible enough to articulate with the adjustable base.

Pressure mapping is a comparative assessment tool. Individual outcomes vary according to body shape, weight, position, mobility and other factors. Mattress selection should form part of an individual assessment where pressure-injury risk is present.

THERAPEDIC® WORLDWIDE

SOY VISCO®

PREMIUM FOAM

EDGE SUPPORT

AUSTRALIAN MADE®

10 years Guarantee

TheraFlex

OT consideration: Encourage adults who have been assessed as being at risk of developing pressure ulcers to change their position frequently and at least every 6 hours. Encourage adults who have been assessed as being at high risk of developing a pressure ulcer to change their position frequently and at least every 4 hours.

T25


Theraflex T25

Feel Rating:
Medium

User Weight:
Max 320kg

Product Highlight:

- Designed for use on adjustable beds
- Flip free design
- Anti slip base layer
- Compatible with waterproof covers



1. Edge Support (65mm Width)

Extra firm edge integrated throughout the mattress to provide complete structural stability and prevent roll-out.

2. Contact Comfort Layer

Soy Visco memory foam, the most natural foam in the market,one of the main ingredients is soy oil from sustainable plantations. It features incredible elasticity and breathability.

3. Premium Support Layer

Engineered with high-density foam, the bottom support layer offers superior structural integrity and weight distribution. It is designed to withstand frequent patient use while minimizing mattress sagging.


4. Airflow Fabric

Promotes breathability keeping the mattress cooler.

Available in all sizes

PRESSURE MAP

Powered & Tested By:
XSENSOR® Technology Corporation
Innovators in Pressure Imaging



TheraFlex T Series Mattress

Standard Mattress

Pressure Gauge

Extreme | High | Medium | Low | Minimum

The way a support surface distributes body weight is an important consideration when selecting a mattress for users with individual positioning or pressure-management requirements.

Theraflex mattress uses layered premium foams designed to distribute load across the mattress surface while remaining flexible enough to articulate with the adjustable base.

Pressure mapping is a comparative assessment tool. Individual outcomes vary according to body shape, weight, position, mobility and other factors. Mattress selection should form part of an individual assessment where pressure-injury risk is present.

THERAPEDIC® WORLDWIDE

SOY VISCO®

PREMIUM FOAM

EDGE SUPPORT

AUSTRALIAN MADE®

10 years Guarantee

TheraFlex

OT consideration: Encourage adults who have been assessed as being at risk of developing pressure ulcers to change their position frequently and at least every 6 hours. Encourage adults who have been assessed as being at high risk of developing a pressure ulcer to change their position frequently and at least every 4 hours.

T30

Available in all sizes


Theraflex T30

Feel Rating:
Plush

User Weight:
Max 320kg


Product Highlight:

- Designed for use on adjustable beds
- Flip free design
- Anti slip base layer
- Compatible with waterproof covers



PRESSURE MAP

Powered & Tested By:
XSENSOR® Technology Corporation
Innovators in Pressure Imaging



TheraFlex T Series Mattress

Standard Mattress

Pressure Gauge

Extreme High Medium Low Minimum

The way a support surface distributes body weight is an important consideration when selecting a mattress for users with individual positioning or pressure-management requirements.

Theraflex mattress uses layered premium foams designed to distribute load across the mattress surface while remaining flexible enough to articulate with the adjustable base.

Pressure mapping is a comparative assessment tool. Individual outcomes vary according to body shape, weight, position, mobility and other factors. Mattress selection should form part of an individual assessment where pressure-injury risk is present.

1. Edge Support (65mm Width)

Extra firm edge integrated throughout the mattress to provide complete structural stability and prevent roll-out.

2. Contact Comfort Layer

CONFIDENCE IN TEXTILES
Tested for harmful substances
according to Oeko-Tex® Standard 100
00000000 Institute

Viscoool®, the phase change material technology provides an open cell structure that guarantees a higher moisture exchange and breathability. It contains 20% soybean oil making it the most natural foam available.

3. Pressure Reduction Layer

Thermal Memory Foam, featuring Graphene - a natural conductor of heat so absorbs body heat to stay cooler and more hygienic.

4. Premium Support Layer

Engineered with high-density foam, the bottom support layer offers superior structural integrity and weight distribution. It is designed to withstand frequent patient use while minimizing mattress sagging.

5. Airflow Fabric

Promotes breathability keeping the mattress cooler.

THERAPEDIC® WORLDWIDE

THERMAL Memory Foam with Graphene

VISCOOOL®

PREMIUM FOAM

EDGE SUPPORT

AUSTRALIAN MADE.

10 years Guarantee

ADVANCED DUAL-TIER PROFILE SPRING TECHNOLOGY FOR ERGONOMIC SPINAL ALIGNMENT

TGA Registered (ARTG: 459568)

Listed on the Australian Register of Therapeutic Goods (TGA Class 1 Medical Device), meeting stringent national safety, therapeutic quality, and manufacturing standards.

Spinal Decompression

Active dual-spring profiling cradles the natural lordotic curve of the lumbar spine, decompressing spinal discs during nocturnal rest.

Clinical Pressure Redistribution

Dynamically diffuses peak pressure points at the shoulders and pelvis to promote healthy vascular circulation and reduce tossing and turning.



TGA MEDICAL DEVICE CLASS 1 (ARTG: 459568)

TheraFlex mattresses are officially registered with the Therapeutic Goods Administration (TGA), the Australian Government Department of Health's regulatory body. Meeting rigorous standards for therapeutic effectiveness, safety, and durability, TheraFlex is designed to support clinical rehabilitation and daily musculoskeletal health.

Therapeutic Indications:

- Postural Correction: Actively assists in maintaining neutral cervical and lumbar alignment in both supine and lateral sleep positions.
- Symptomatic Relief: Engineered to alleviate tension in the lower back (lumbar), reduce sciatica discomfort, and ease morning stiffness.
- Vascular Support: Prevents localized capillary compression, helping to eliminate numbness and tingling in high-load areas.clinical rehabilitation and daily musculoskeletal health.

DUAL-TIER PROFILE SPRING SYSTEM

- Targeted Anatomical Contouring:
  Engineered with varying spring tensions to conform precisely to human body biomechanics, providing gentle immersion for broader shoulder and pelvic bones while delivering reinforced upward thrust to the lumbar arch.
- Gradual Multi-Layer Deflection:
  The top soft-gauge coils absorb initial kinetic impact for immediate surface pressure relief; the high-tensile base layer acts as an ergonomic stabilizer to prevent spine sag.
- Independent Motion Isolation:
  Individually pocketed coils respond to micro-movements independently, ensuring zero cross-mattress disturbance and uninterrupted deep-sleep cycles.

Scan to Verify:

Verify TGA Medical Device Registration (ARTG ID: 459568) via the official Australian Government Therapeutic Goods Administration portal.

TheraFlex

OT consideration: Encourage adults who have been assessed as being at risk of developing pressure ulcers to change their position frequently and at least every 6 hours. Encourage adults who have been assessed as being at high risk of developing a pressure ulcer to change their position frequently and at least every 4 hours.

PROFILE THERAPY COMFORT

Available in all sizes

Profile Therapy Comfort

Feel Rating:
Medium


User Weight:

Max 320kg

Total Height:
27cm

Product Highlight:

- Designed for use on adjustable beds
- Flip free design
- Anti slip base layer
- Compatible with waterproof covers



1. Adaptive Fabric

The most dynamic and innovative bedding material available.Tests performed independently by North Carolina State University and RMIT University in Melbourne proves Adaptive fabric improves the wicking and evaporation of humidity away from the body and regulates the temperature in a sleep environment.

2. Premium Foam

High quality premium comfort layer provides extra comfort. High breathable and medium feel.

3. Profile Spring

The top springs are soft and bottom springs firm giving incredible gradual pressure relief. It's designed to profile to your body shape maximising support where you need it.

4. Hinged Probox


Superior 90mm wide side support around the perimeter of the mattress. Provides edge to edge sleeping surface without roll out. Hinged walls maximize the flexibility.

TGA MEDICAL DEVICE CLASS 1 (ARTG: 459568)

Designed and engineered to contour to the natural spinal curvature and provide extra support in lower back area

Unique 2 layered pocket spring system has the top layer softer than the bottom layer for gradual pressure relief and curved for spinal alignment

Pressure Map



TheraFlex T Series Mattress | Standard Mattress

Pressure Gauge

Extreme | Low | High | Minimum | Medium

The way a support surface distributes body weight is an important consideration when selecting a mattress for users with individual positioning or pressure-management requirements.

Theraflex mattress uses layered premium foams designed to distribute load across the mattress surface while remaining flexible enough to articulate with the adjustable base.

Pressure mapping is a comparative assessment tool. Individual outcomes vary according to body shape, weight, position, mobility and other factors. Mattress selection should form part of an individual assessment where pressure-injury risk is present.

Powered & Tested By:
XSENSOR® Technology Corporation
Innovators in Pressure Imaging

THERAPEDIC® WORLDWIDE

PREMIUM FOAM | PROFILE SPRING | HINGED PROBOX® Edge Support System | AUSTRALIAN MADE | 10 years Guarantee

TheraFlex

OT consideration: Encourage adults who have been assessed as being at risk of developing pressure ulcers to change their position frequently and at least every 6 hours. Encourage adults who have been assessed as being at high risk of developing a pressure ulcer to change their position frequently and at least every 4 hours.

PROFILE THERAPY LUXURY

Available in all sizes

Profile Therapy Luxury

Feel Rating:
Plush


User Weight:

Max 320kg

Total Height:
29cm

Product Highlight:

- Designed for use on adjustable beds
- Flip free design
- Anti slip base layer
- Compatible with waterproof covers



1. Adaptive Fabric

The most dynamic and innovative bedding material available.Tests performed independently by North Carolina State University and RMIT University in Melbourne proves Adaptive fabric improves the wicking and evaporation of humidity away from the body and regulates the temperature in a sleep environment.

2. Cloudflex

tested by RMIT UNIVERSITY

Latest foam technology backed by RMIT research study into quality performance- proved to have superior pressure redistribution and heat dissipation.

3. Premium Foam

High quality premium comfort layer provides extra comfort. High breathable and plush feel.

5. Profile Spring

The top springs are soft and bottom springs firm giving incredible gradual pressure relief. It's designed to profile to your body shape maximising support where you need it.

4. Hinged Probox


Superior 90mm wide side support around the perimeter of the mattress. Provides edge to edge sleeping surface without roll out. Hinged walls maximize the flexibility.

TGA MEDICAL DEVICE CLASS 1 (ARTG: 459568)

Designed and engineered to contour to the natural spinal curvature and provide extra support in lower back area

Unique 2 layered pocket spring system has the top layer softer than the bottom layer for gradual pressure relief and curved for spinal alignment

Pressure Map



TheraFlex T Series Mattress | Standard Mattress

Pressure Gauge

Extreme | Low

High | Minimum

Medium

The way a support surface distributes body weight is an important consideration when selecting a mattress for users with individual positioning or pressure-management requirements.

Theraflex mattress uses layered premium foams designed to distribute load across the mattress surface while remaining flexible enough to articulate with the adjustable base.

Pressure mapping is a comparative assessment tool. Individual outcomes vary according to body shape, weight, position, mobility and other factors. Mattress selection should form part of an individual assessment where pressure-injury risk is present.

Powered & Tested By:
XSENSOR® Technology Corporation
Innovators in Pressure Imaging

THERAPEDIC® WORLDWIDE

Cloudflex OEKO-TEX® STANDARD 100 APPROVED

PREMIUM FOAM

PROFILE SPRING

HINGED PROBOX® Edge Support System

AUSTRALIAN MADE

10 years Guarantee